

TNFR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00777-2

Capture Antibody Information

Catalog Number:
83101-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230328E8
Reactivity:
human
GenBank:
BC052977

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 1B
Gene ID:
7133

Detection Antibody Information

Catalog Number:
98010-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230328B3
Reactivity:
human
GenBank:
BC052977

Conjugate:
Unconjugated
Full name:
tumor necrosis factor receptor superfamily, member 1B
Gene ID:
7133

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00777-2 targets TNFR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFR2 Recombinant antibody, PBS Only (Capture) 83101-3-PBS (230328E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: Anti-Human TNFR2/CD120b Rabbit Recombinant Antibody, PBS Only 98010-1-PBS (230328B3). 100 µg. Concentration 1 mg/ml.

Alternative TNFR2 matched antibody pairs: MP00777-1, MP00777-3, MP00777-4, MP50139-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00777-2, TNFRSF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83101-3-PBS. Detection antibody: 98010-1-PBS. Standard: Eg0086. Range: 0.625-20 ng/mL.