

# TNFR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00777-3

## Capture Antibody Information

**Catalog Number:**  
83101-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230328F5  
**Reactivity:**  
human  
**GenBank:**  
BC052977

**Conjugate:**  
Unconjugated  
**Full name:**  
tumor necrosis factor receptor superfamily, member 1B  
**Gene ID:**  
7133

## Detection Antibody Information

**Catalog Number:**  
98010-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
230328B3  
**Reactivity:**  
human  
**GenBank:**  
BC052977

**Conjugate:**  
Unconjugated  
**Full name:**  
tumor necrosis factor receptor superfamily, member 1B  
**Gene ID:**  
7133

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00777-3 targets TNFR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFR2 Recombinant antibody, PBS Only (Capture) 83101-4-PBS (230328F5). 100 µg. Concentration 1 mg/ml.

Detection antibody: Anti-Human TNFR2/CD120b Rabbit Recombinant Antibody, PBS Only 98010-1-PBS (230328B3). 100 µg. Concentration 1 mg/ml.

Alternative TNFR2 matched antibody pairs: MP00777-1, MP00777-2, MP00777-4, MP50139-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00777-3, TNFRSF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83101-4-PBS. Detection antibody: 98010-1-PBS. Standard: Eg0086. Range: 0.625-20 ng/mL.