

Aquaporin 4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00781-2

Capture Antibody Information

Catalog Number:
82296-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240856B2

Reactivity:
human

GenBank:
BC022286

Immunogen Catalog Number:
Ag9561

Conjugate:
Unconjugated

Full name:
Aquaporin 4

Gene ID:
361

Detection Antibody Information

Catalog Number:
82296-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240856A9

Reactivity:
human

GenBank:
BC022286

Immunogen Catalog Number:
Ag9561

Conjugate:
Unconjugated

Full name:
Aquaporin 4

Gene ID:
361

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00781-2 targets Aquaporin 4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Aquaporin 4 Recombinant antibody, PBS Only (Capture/Detector) 82296-3-PBS (240856B2). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Aquaporin 4 Recombinant antibody, PBS Only (Detector) 82296-2-PBS (240856A9). 100 μ g. Concentration 1 mg/ml.

Alternative Aquaporin 4 matched antibody pairs: MP00781-1, MP00781-3, MP00781-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00781-2, AQP4/Aquaporin 4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82296-3-PBS. Detection antibody: 82296-2-PBS. Standard: Ag9561. Range: 0.313-40 ng/mL