

FRMD3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00787-1

Capture Antibody Information

Catalog Number: 83850-3-PBS	Clone ID: 240967G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FERM domain containing 3
Isotype: IgG	GenBank: BC023560	Gene ID: 257019
Purification Method: Protein A purification	Immunogen Catalog Number: Ag13235	

Detection Antibody Information

Catalog Number: 83850-2-PBS	Clone ID: 240967D9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: FERM domain containing 3
Isotype: IgG	GenBank: BC023560	Gene ID: 257019
Purification Method: Protein A purification	Immunogen Catalog Number: Ag13235	

Applications

Tested Applications: Cytometric bead array	Range: 0.78-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00787-1 targets FRMD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FRMD3 Recombinant antibody, PBS Only (Capture) 83850-3-PBS (240967G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: FRMD3 Recombinant antibody, PBS Only (Detector) 83850-2-PBS (240967D9). 100 µg. Concentration 1 mg/mL.

Alternative FRMD3 matched antibody pairs: MP00787-2, MP00787-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00787-1, FRMD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83850-3-PBS. Detection antibody: 83850-2-PBS. Standard: Ag13235. Range: 0.78-100 ng/mL