

CLCA1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00794-2

Capture Antibody Information

Catalog Number:
83829-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240941C4
Reactivity:
human
GenBank:
BC156805
Immunogen Catalog Number:
Ag18176

Conjugate:
Unconjugated
Full name:
chloride channel accessory 1
Gene ID:
1179

Detection Antibody Information

Catalog Number:
83829-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240941G6
Reactivity:
human
GenBank:
BC156805
Immunogen Catalog Number:
Ag18176

Conjugate:
Unconjugated
Full name:
chloride channel accessory 1
Gene ID:
1179

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00794-2 targets CLCA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CLCA1 Recombinant antibody, PBS Only (Capture) 83829-2-PBS (240941C4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CLCA1 Recombinant antibody, PBS Only (Detector) 83829-5-PBS (240941G6). 100 µg. Concentration 1 mg/ml.

Alternative CLCA1 matched antibody pairs: MP00794-1, MP00794-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00794-2, CLCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83829-2-PBS. Detection antibody: 83829-5-PBS. Standard: Ag18176. Range: 0.625-80 ng/mL