

DLL3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00798-1

Capture Antibody Information

Catalog Number:
83834-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240956D5

Reactivity:
human

GenBank:
BC000218

Immunogen Catalog Number:
Ag0579

Conjugate:
Unconjugated

Full name:
delta-like 3 (Drosophila)

Gene ID:
10683

Detection Antibody Information

Catalog Number:
83834-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240956D3

Reactivity:
human

GenBank:
BC000218

Immunogen Catalog Number:
Ag0579

Conjugate:
Unconjugated

Full name:
delta-like 3 (Drosophila)

Gene ID:
10683

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00798-1 targets DLL3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DLL3 Recombinant antibody, PBS Only (Capture) 83834-2-PBS (240956D5). 100 µg. Concentration 1 mg/ml.

Detection antibody: DLL3 Recombinant antibody, PBS Only (Detector) 83834-1-PBS (240956D3). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00798-1, DLL3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83834-2-PBS. Detection antibody: 83834-1-PBS. Standard: Ag0579. Range: 0.313-40 ng/mL