

For Research Use Only

# GADD45GIP1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00800-1

## Capture Antibody Information

Catalog Number:  
83823-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240947A9

Reactivity:  
human

GenBank:  
BC013027

Immunogen Catalog Number:  
Ag9299

Conjugate:  
Unconjugated

Full name:  
growth arrest and DNA-damage-inducible, gamma interacting protein 1

Gene ID:  
90480

## Detection Antibody Information

Catalog Number:  
83823-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240947D12

Reactivity:  
human

GenBank:  
BC013027

Immunogen Catalog Number:  
Ag9299

Conjugate:  
Unconjugated

Full name:  
growth arrest and DNA-damage-inducible, gamma interacting protein 1

Gene ID:  
90480

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00800-1 targets GADD45GIP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GADD45GIP1 Recombinant antibody, PBS Only (Capture) 83823-3-PBS (240947A9). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: GADD45GIP1 Recombinant antibody, PBS Only (Detector) 83823-4-PBS (240947D12). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative GADD45GIP1 matched antibody pairs: MP00800-2, MP00800-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00800-1, GADD45GIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83823-3-PBS. Detection antibody: 83823-4-PBS. Standard: Ag9299. Range: 0.78-100 ng/mL.