

FBN1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00807-2

Capture Antibody Information

Catalog Number:
83836-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240960G1

Reactivity:
human

GenBank:
BC146854

Immunogen Catalog Number:
Ag30283

Conjugate:
Unconjugated

Full name:
fibrillin 1

Gene ID:
2200

Detection Antibody Information

Catalog Number:
83836-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240960E5

Reactivity:
human

GenBank:
BC146854

Immunogen Catalog Number:
Ag30283

Conjugate:
Unconjugated

Full name:
fibrillin 1

Gene ID:
2200

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00807-2 targets FBN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FBN1 Recombinant antibody, PBS Only (Capture) 83836-3-PBS (240960G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: FBN1 Recombinant antibody, PBS Only (Detector) 83836-1-PBS (240960E5). 100 µg. Concentration 1 mg/ml.

Alternative FBN1 matched antibody pairs: MP00807-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00807-2, Fibrillin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83836-3-PBS. Detection antibody: 83836-1-PBS. Standard: Ag30283. Range: 0.313-40 ng/mL.