

For Research Use Only

OLIG2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00809-1

Capture Antibody Information

Catalog Number:
82806-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240839D6
Reactivity:
human
GenBank:
BC047511
Immunogen Catalog Number:
Ag5089

Conjugate:
Unconjugated
Full name:
oligodendrocyte lineage transcription factor 2
Gene ID:
10215

Detection Antibody Information

Catalog Number:
82806-6-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240839F5
Reactivity:
human
GenBank:
BC047511
Immunogen Catalog Number:
Ag5089

Conjugate:
Unconjugated
Full name:
oligodendrocyte lineage transcription factor 2
Gene ID:
10215

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00809-1 targets OLIG2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OLIG2 Recombinant antibody, PBS Only (Capture) 82806-5-PBS (240839D6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: OLIG2 Recombinant antibody, PBS Only (Detector) 82806-6-PBS (240839F5). 100 μ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00809-1, OLIG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82806-5-PBS. Detection antibody: 82806-6-PBS. Standard: Ag5089. Range: 0.625-80 ng/mL