

For Research Use Only

VWF Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00810-1

Capture Antibody Information

Catalog Number:
83854-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240867H2
Reactivity:
human
GenBank:
Immunogen Catalog Number:
Ag25578

Conjugate:
Unconjugated
Full name:
von Willebrand factor
Gene ID:
7450

Detection Antibody Information

Catalog Number:
83854-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240867D7
Reactivity:
human
GenBank:
Immunogen Catalog Number:
Ag25578

Conjugate:
Unconjugated
Full name:
von Willebrand factor
Gene ID:
7450

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00810-1 targets VWF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: vwf Recombinant antibody, PBS Only (Capture) 83854-3-PBS (240867H2). 100 μ g. Concentration 1 mg/mL.
Detection antibody: vwf Recombinant antibody, PBS Only (Detector) 83854-1-PBS (240867D7). 100 μ g. Concentration 1 mg/mL.
Alternative VWF matched antibody pairs: MP00810-2, MP00810-3
Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.
Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00810-1, vwf Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83854-3-PBS. Detection antibody: 83854-1-PBS. Standard: Ag25578. Range: 0.313-40 ng/mL