

TFPI Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00812-1

Capture Antibody Information

Catalog Number:
83838-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240977C12
Reactivity:
human
GenBank:
NM_006287.6

Conjugate:
Unconjugated
Full name:
tissue factor pathway inhibitor (lipoprotein-associated coagulation inhibitor)
Gene ID:
7035

Detection Antibody Information

Catalog Number:
83838-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240977D6
Reactivity:
human
GenBank:
NM_006287.6

Conjugate:
Unconjugated
Full name:
tissue factor pathway inhibitor (lipoprotein-associated coagulation inhibitor)
Gene ID:
7035

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00812-1 targets TFPI in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TFPI Recombinant antibody, PBS Only (Capture) 83838-1-PBS (240977C12). 100 μ g. Concentration 1 mg/ml.

Detection antibody: TFPI Recombinant antibody, PBS Only (Detector) 83838-3-PBS (240977D6). 100 μ g. Concentration 1 mg/ml.

Alternative TFPI matched antibody pairs: MP00812-2, MP00812-3, MP00812-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00812-1, TFPI Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83838-1-PBS. Detection antibody: 83838-3-PBS. Standard: Eg0855. Range: 0.781-100 ng/mL.