

TRMT6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00816-2

Capture Antibody Information

Catalog Number: 83843-1-PBS	Clone ID: 240911A1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tRNA methyltransferase 6 homolog (S. cerevisiae)
Isotype: IgG	GenBank: BC001262	Gene ID: 51605
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10176	

Detection Antibody Information

Catalog Number: 83843-4-PBS	Clone ID: 240911F12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: tRNA methyltransferase 6 homolog (S. cerevisiae)
Isotype: IgG	GenBank: BC001262	Gene ID: 51605
Purification Method: Protein A purification	Immunogen Catalog Number: Ag10176	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00816-2 targets TRMT6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRMT6 Recombinant antibody, PBS Only (Capture) 83843-1-PBS (240911A1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRMT6 Recombinant antibody, PBS Only (Capture/Detector) 83843-4-PBS (240911F12). 100 µg. Concentration 1 mg/ml.

Alternative TRMT6 matched antibody pairs: MP00816-1, MP00816-3


Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00816-2, TRMT6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83843-1-PBS. Detection antibody: 83843-4-PBS. Standard: Ag10176. Range: 0.625-80 ng/mL