

EXOSC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00826-3

Capture Antibody Information

Catalog Number: 83875-4-PBS	Clone ID: 240938F3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: exosome component 1
Isotype: IgG	GenBank: BC022067	Gene ID: 51013
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3286	

Detection Antibody Information

Catalog Number: 83875-2-PBS	Clone ID: 240938E9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: exosome component 1
Isotype: IgG	GenBank: BC022067	Gene ID: 51013
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3286	

Applications

Tested Applications: Cytometric bead array	Range: 0.78-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00826-3 targets EXOSC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EXOSC1 Recombinant antibody, PBS Only (Capture) 83875-4-PBS (240938F3). 100 µg. Concentration 1 mg/ml.

Detection antibody: EXOSC1 Recombinant antibody, PBS Only (Capture/Detector) 83875-2-PBS (240938E9). 100 µg. Concentration 1 mg/ml.

Alternative EXOSC1 matched antibody pairs: MP00826-1, MP00826-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00826-3, EXOSC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83875-4-PBS. Detection antibody: 83875-2-PBS. Standard: Ag3286. Range: 0.78-100 ng/mL.