

For Research Use Only

NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00833-2

Capture Antibody Information

Catalog Number:
83882-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240809E8
Reactivity:
human
GenBank:
BC052582

Conjugate:
Unconjugated
Full name:
natural cytotoxicity triggering receptor 3
Gene ID:
259197

Detection Antibody Information

Catalog Number:
83882-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240809C7
Reactivity:
human
GenBank:
BC052582

Conjugate:
Unconjugated
Full name:
natural cytotoxicity triggering receptor 3
Gene ID:
259197

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00833-2 targets NCR3/NKp30 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Capture) 83882-3-PBS (240809E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Detector) 83882-2-PBS (240809C7). 100 µg. Concentration 1 mg/ml.

Alternative NCR3/NKp30 matched antibody pairs: MP00833-1, MP00833-3, MP00833-4, MP50119-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00833-2, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-3-PBS. Detection antibody: 83882-2-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.