

# NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00833-3

## Capture Antibody Information

**Catalog Number:**  
83882-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240809E8  
**Reactivity:**  
human  
**GenBank:**  
BC052582

**Conjugate:**  
Unconjugated  
**Full name:**  
natural cytotoxicity triggering receptor 3  
**Gene ID:**  
259197

## Detection Antibody Information

**Catalog Number:**  
83882-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240809A9  
**Reactivity:**  
human  
**GenBank:**  
BC052582

**Conjugate:**  
Unconjugated  
**Full name:**  
natural cytotoxicity triggering receptor 3  
**Gene ID:**  
259197

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00833-3 targets NCR3/NKp30 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Capture) 83882-3-PBS (240809E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: NCR3/NKp30 Recombinant antibody, PBS Only (Detector) 83882-1-PBS (240809A9). 100 µg. Concentration 1 mg/ml.

Alternative NCR3/NKp30 matched antibody pairs: MP00833-1, MP00833-2, MP00833-4, MP50119-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00833-3, NCR3/NKp30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83882-3-PBS. Detection antibody: 83882-1-PBS. Standard: Eg31663. Range: 0.313-40 ng/mL.