

For Research Use Only

SIK1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00838-1

Capture Antibody Information

Catalog Number: 83897-3-PBS	Clone ID: 240961F9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: salt-inducible kinase 1
Isotype: IgG	GenBank: BC038504	Gene ID: 150094
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28481	

Detection Antibody Information

Catalog Number: 83897-1-PBS	Clone ID: 240961C2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: salt-inducible kinase 1
Isotype: IgG	GenBank: BC038504	Gene ID: 150094
Purification Method: Protein A purification	Immunogen Catalog Number: Ag28481	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00838-1 targets SIK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SIK1 Recombinant antibody, PBS Only (Capture) 83897-3-PBS (240961F9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SIK1 Recombinant antibody, PBS Only (Detector) 83897-1-PBS (240961C2). 100 μ g. Concentration 1 mg/ml.

Alternative SIK1 matched antibody pairs: MP00838-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00838-1, SIK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83897-3-PBS. Detection antibody: 83897-1-PBS. Standard: Ag28481. Range: 0.313-40 ng/mL