

VISTA Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00865-2

Capture Antibody Information

Catalog Number:
83919-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241003F4

Reactivity:
human

GenBank:
BC020568

Conjugate:
Unconjugated

Full name:
chromosome 10 open reading frame 54

Gene ID:
64115

Detection Antibody Information

Catalog Number:
83919-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241003A12

Reactivity:
human

GenBank:
BC020568

Conjugate:
Unconjugated

Full name:
chromosome 10 open reading frame 54

Gene ID:
64115

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00865-2 targets VISTA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: VISTA Recombinant antibody, PBS Only (Capture) 83919-3-PBS (241003F4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: VISTA Recombinant antibody, PBS Only (Detector) 83919-2-PBS (241003A12). 100 μ g. Concentration 1 mg/ml.

Alternative VISTA matched antibody pairs: MP00865-1, MP50058-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00865-2, VISTA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83919-3-PBS. Detection antibody: 83919-2-PBS. Standard: Eg0017. Range: 0.156-20 ng/mL