

For Research Use Only

Mouse IL23R Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00866-1

Capture Antibody Information

Catalog Number:
83918-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240989B3
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
interleukin 23 receptor
Gene ID:
209590

Detection Antibody Information

Catalog Number:
83918-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240989G1
Reactivity:
mouse
GenBank:

Conjugate:
Unconjugated
Full name:
interleukin 23 receptor
Gene ID:
209590

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00866-1 targets IL23R in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse IL23R Recombinant antibody, PBS Only (Capture) 83918-1-PBS (240989B3). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse IL23R Recombinant antibody, PBS Only (Detector) 83918-2-PBS (240989G1). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00866-1, MOUSE IL23R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83918-1-PBS. Detection antibody: 83918-2-PBS. Standard: Eg0808. Range: 0.625-40 ng/mL.