

PHF19 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00868-2

Capture Antibody Information

Catalog Number: 83892-5-PBS	Clone ID: 241047G10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: PHD finger protein 19
Isotype: IgG	GenBank: BC022374	Gene ID: 26147
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2511	

Detection Antibody Information

Catalog Number: 83892-3-PBS	Clone ID: 241047D7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: PHD finger protein 19
Isotype: IgG	GenBank: BC022374	Gene ID: 26147
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2511	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00868-2 targets PHF19 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PHF19 Recombinant antibody, PBS Only (Capture) 83892-5-PBS (241047G10). 100 µg. Concentration 1 mg/ml.

Detection antibody: PHF19 Recombinant antibody, PBS Only (Detector) 83892-3-PBS (241047D7). 100 µg. Concentration 1 mg/ml.

Alternative PHF19 matched antibody pairs: MP00868-1, MP00868-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00868-2, PHF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83892-5-PBS. Detection antibody: 83892-3-PBS. Standard: Ag2511. Range: 0.313-40 ng/mL