

Securin Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00876-3

Capture Antibody Information

Catalog Number:
83981-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241090G10
Reactivity:
human
GenBank:
BC101834
Immunogen Catalog Number:
Ag12670

Conjugate:
Unconjugated
Full name:
pituitary tumor-transforming 1
Gene ID:
9232

Detection Antibody Information

Catalog Number:
83981-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241090F9
Reactivity:
human
GenBank:
BC101834
Immunogen Catalog Number:
Ag12670

Conjugate:
Unconjugated
Full name:
pituitary tumor-transforming 1
Gene ID:
9232

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00876-3 targets Securin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Securin Recombinant antibody, PBS Only (Capture/Detector) 83981-4-PBS (241090G10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Securin Recombinant antibody, PBS Only (Detector) 83981-3-PBS (241090F9). 100 μ g. Concentration 1 mg/mL.

Alternative Securin matched antibody pairs: MP00876-1, MP00876-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00876-3, Securin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83981-4-PBS. Detection antibody: 83981-3-PBS. Standard: Ag12670. Range: 1.56-100 ng/mL