

TPRG1L Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00883-1

Capture Antibody Information

Catalog Number:
83989-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241070A1
Reactivity:
human
GenBank:
BC019034
Immunogen Catalog Number:
Ag22585

Conjugate:
Unconjugated
Full name:
tumor protein p63 regulated 1-like
Gene ID:
127262

Detection Antibody Information

Catalog Number:
83989-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241070B8
Reactivity:
human
GenBank:
BC019034
Immunogen Catalog Number:
Ag22585

Conjugate:
Unconjugated
Full name:
tumor protein p63 regulated 1-like
Gene ID:
127262

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00883-1 targets TPRG1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TPRG1L Recombinant antibody, PBS Only (Capture) 83989-1-PBS (241070A1). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TPRG1L Recombinant antibody, PBS Only (Detector) 83989-3-PBS (241070B8). 100 μ g. Concentration 1 mg/mL.

Alternative TPRG1L matched antibody pairs: MP00883-2, MP00883-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00883-1, TPRG1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83989-1-PBS. Detection antibody: 83989-3-PBS. Standard: Ag22585. Range: 0.625-80 ng/mL