

For Research Use Only

AGXT2L1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00885-3

Capture Antibody Information

Catalog Number:
83984-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241119G4
Reactivity:
human
GenBank:
BC022526
Immunogen Catalog Number:
Ag10584

Conjugate:
Unconjugated
Full name:
alanine-glyoxylate aminotransferase 2-like 1
Gene ID:
64850

Detection Antibody Information

Catalog Number:
83984-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241119D3
Reactivity:
human
GenBank:
BC022526
Immunogen Catalog Number:
Ag10584

Conjugate:
Unconjugated
Full name:
alanine-glyoxylate aminotransferase 2-like 1
Gene ID:
64850

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00885-3 targets AGXT2L1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AGXT2L1 Recombinant antibody, PBS Only (Capture/Detector) 83984-4-PBS (241119G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: AGXT2L1 Recombinant antibody, PBS Only (Detector) 83984-3-PBS (241119D3). 100 µg. Concentration 1 mg/ml.

Alternative AGXT2L1 matched antibody pairs: MP00885-1, MP00885-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00885-3, AGXT2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83984-4-PBS. Detection antibody: 83984-3-PBS. Standard: Ag10584. Range: 0.625-80 ng/mL.