

CRMP2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00886-1

Capture Antibody Information

Catalog Number:
83990-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241088A5
Reactivity:
human
GenBank:
BC067109
Immunogen Catalog Number:
Ag5987

Conjugate:
Unconjugated
Full name:
dihydropyrimidinase-like 2
Gene ID:
1808

Detection Antibody Information

Catalog Number:
83990-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241088C6
Reactivity:
human
GenBank:
BC067109
Immunogen Catalog Number:
Ag5987

Conjugate:
Unconjugated
Full name:
dihydropyrimidinase-like 2
Gene ID:
1808

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00886-1 targets CRMP2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CRMP2 Recombinant antibody, PBS Only (Capture) 83990-1-PBS (241088A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: CRMP2 Recombinant antibody, PBS Only (Detector) 83990-2-PBS (241088C6). 100 µg. Concentration 1 mg/ml.

Alternative CRMP2 matched antibody pairs: MP00886-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00886-1, DPYSL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83990-1-PBS. Detection antibody: 83990-2-PBS. Standard: Ag5987. Range: 6.25-100 ng/mL