

PI16 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00888-2

Capture Antibody Information

Catalog Number: 83985-1-PBS	Clone ID: 241014B5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: peptidase inhibitor 16
Isotype: IgG	GenBank: BC022399	Gene ID: 221476
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2925	

Detection Antibody Information

Catalog Number: 83985-2-PBS	Clone ID: 241014B11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: peptidase inhibitor 16
Isotype: IgG	GenBank: BC022399	Gene ID: 221476
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2925	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00888-2 targets PI16 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PI16 Recombinant antibody, PBS Only (Capture) 83985-1-PBS (241014B5). 100 µg. Concentration 1 mg/mL.

Detection antibody: PI16 Recombinant antibody, PBS Only (Detector) 83985-2-PBS (241014B11). 100 µg. Concentration 1 mg/mL.

Alternative PI16 matched antibody pairs: MP00888-1, MP00888-3, MP00888-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00888-2, PI16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83985-1-PBS. Detection antibody: 83985-2-PBS. Standard: Ag2925. Range: 0.313-40 ng/mL