

For Research Use Only

# SNAT2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00889-1

## Capture Antibody Information

Catalog Number:  
83968-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241059A10

Reactivity:  
human

GenBank:  
BC040342

Immunogen Catalog Number:  
Ag23082

Conjugate:  
Unconjugated

Full name:  
solute carrier family 38, member 2

Gene ID:  
54407

## Detection Antibody Information

Catalog Number:  
83968-2-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241059D8

Reactivity:  
human

GenBank:  
BC040342

Immunogen Catalog Number:  
Ag23082

Conjugate:  
Unconjugated

Full name:  
solute carrier family 38, member 2

Gene ID:  
54407

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00889-1 targets SNAT2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC38A2 Recombinant antibody, PBS Only (Capture/Detector) 83968-1-PBS (241059A10). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC38A2 Recombinant antibody, PBS Only (Detector) 83968-2-PBS (241059D8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00889-1, SLC38A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83968-1-PBS. Detection antibody: 83968-2-PBS. Standard: Ag23082. Range: 0.313-40 ng/mL.