

For Research Use Only

QSOX1/Quiescin Q6 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00903-1

Capture Antibody Information

Catalog Number:
83952-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240946D5
Reactivity:
human
GenBank:
BC017692
Immunogen Catalog Number:
Ag3396

Conjugate:
Unconjugated
Full name:
quiescin Q6 sulfhydryl oxidase 1
Gene ID:
5768

Detection Antibody Information

Catalog Number:
83952-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240946B3
Reactivity:
human
GenBank:
BC017692
Immunogen Catalog Number:
Ag3396

Conjugate:
Unconjugated
Full name:
quiescin Q6 sulfhydryl oxidase 1
Gene ID:
5768

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00903-1 targets QSOX1/Quiescin Q6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: QSOX1/Quiescin Q6 Recombinant antibody, PBS Only (Capture) 83952-2-PBS (240946D5). 100 µg. Concentration 1 mg/ml.

Detection antibody: QSOX1/Quiescin Q6 Recombinant antibody, PBS Only (Detector) 83952-1-PBS (240946B3). 100 µg. Concentration 1 mg/ml.

Alternative QSOX1/Quiescin Q6 matched antibody pairs: MP00903-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00903-1, QSOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83952-2-PBS. Detection antibody: 83952-1-PBS. Standard: Ag3396. Range: 1.563-100 ng/mL.