

For Research Use Only

C10orf33 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00905-1

Capture Antibody Information

Catalog Number:
83954-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240984F9
Reactivity:
human
GenBank:
BC006131
Immunogen Catalog Number:
Ag14701

Conjugate:
Unconjugated
Full name:
chromosome 10 open reading frame 33
Gene ID:
84795

Detection Antibody Information

Catalog Number:
83954-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240984A12
Reactivity:
human
GenBank:
BC006131
Immunogen Catalog Number:
Ag14701

Conjugate:
Unconjugated
Full name:
chromosome 10 open reading frame 33
Gene ID:
84795

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00905-1 targets C10orf33 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C10orf33 Recombinant antibody, PBS Only (Capture/Detector) 83954-2-PBS (240984F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: C10orf33 Recombinant antibody, PBS Only (Detector) 83954-1-PBS (240984A12). 100 µg. Concentration 1 mg/ml.

Alternative C10orf33 matched antibody pairs: MP00905-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00905-1, C10orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83954-2-PBS. Detection antibody: 83954-1-PBS. Standard: Ag14701. Range: 1.25-80 ng/mL.