

For Research Use Only

# TLE3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00909-1

## Capture Antibody Information

**Catalog Number:**  
83955-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240680E4  
**Reactivity:**  
human  
**GenBank:**  
BC015729  
**Immunogen Catalog Number:**  
Ag1930

**Conjugate:**  
Unconjugated  
**Full name:**  
transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)  
**Gene ID:**  
7090

## Detection Antibody Information

**Catalog Number:**  
83955-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240680A9  
**Reactivity:**  
human  
**GenBank:**  
BC015729  
**Immunogen Catalog Number:**  
Ag1930

**Conjugate:**  
Unconjugated  
**Full name:**  
transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)  
**Gene ID:**  
7090

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00909-1 targets TLE3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TLE3 Recombinant antibody, PBS Only (Capture) 83955-3-PBS (240680E4). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: TLE3 Recombinant antibody, PBS Only (Detector) 83955-1-PBS (240680A9). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative TLE3 matched antibody pairs: MP00909-2, MP00909-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply. Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00909-1, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-3-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL