

TOB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00930-1

Capture Antibody Information

Catalog Number: 83944-4-PBS	Clone ID: 241043G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: transducer of ERBB2, 1
Isotype: IgG	GenBank: BC031406	Gene ID: 10140
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6714	

Detection Antibody Information

Catalog Number: 83944-3-PBS	Clone ID: 241043C11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: transducer of ERBB2, 1
Isotype: IgG	GenBank: BC031406	Gene ID: 10140
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6714	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00930-1 targets TOB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TOB1 Recombinant antibody, PBS Only (Capture) 83944-4-PBS (241043G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TOB1 Recombinant antibody, PBS Only (Detector) 83944-3-PBS (241043C11). 100 µg. Concentration 1 mg/ml.

Alternative TOB1 matched antibody pairs: MP00930-2, MP00930-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00930-1, TOB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83944-4-PBS. Detection antibody: 83944-3-PBS. Standard: Ag6714. Range: 0.156-20 ng/mL