

For Research Use Only

# TOB1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00930-3

## Capture Antibody Information

<b>Catalog Number:</b> 83944-1-PBS	<b>Clone ID:</b> 241043A5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> transducer of ERBB2, 1
<b>Isotype:</b> IgG	<b>GenBank:</b> BC031406	<b>Gene ID:</b> 10140
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag6714	

## Detection Antibody Information

<b>Catalog Number:</b> 83944-2-PBS	<b>Clone ID:</b> 241043A9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> transducer of ERBB2, 1
<b>Isotype:</b> IgG	<b>GenBank:</b> BC031406	<b>Gene ID:</b> 10140
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag6714	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.313-20 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00930-3 targets TOB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TOB1 Recombinant antibody, PBS Only (Capture) 83944-1-PBS (241043A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: TOB1 Recombinant antibody, PBS Only (Detector) 83944-2-PBS (241043A9). 100 µg. Concentration 1 mg/ml.

Alternative TOB1 matched antibody pairs: MP00930-1, MP00930-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00930-3, TOB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83944-1-PBS. Detection antibody: 83944-2-PBS. Standard: Ag6714. Range: 0.313-20 ng/mL