

ACSF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00952-3

Capture Antibody Information

Catalog Number:
84033-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241208C1
Reactivity:
human
GenBank:
BC012053
Immunogen Catalog Number:
Ag9199

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase family member 2
Gene ID:
80221

Detection Antibody Information

Catalog Number:
84033-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241208D9
Reactivity:
human
GenBank:
BC012053
Immunogen Catalog Number:
Ag9199

Conjugate:
Unconjugated
Full name:
acyl-CoA synthetase family member 2
Gene ID:
80221

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00952-3 targets ACSF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ACSF2 Recombinant antibody, PBS Only (Capture) 84033-1-PBS (241208C1). 100 µg. Concentration 1 mg/ml.

Detection antibody: ACSF2 Recombinant antibody, PBS Only (Detector) 84033-3-PBS (241208D9). 100 µg. Concentration 1 mg/ml.

Alternative ACSF2 matched antibody pairs: MP00952-1, MP00952-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00952-3, ACSF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84033-1-PBS. Detection antibody: 84033-3-PBS. Standard: Ag9199. Range: 2.5-80 ng/mL.