

ABR Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00961-1

Capture Antibody Information

Catalog Number: 84011-1-PBS	Clone ID: 241121B4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: active BCR-related gene
Isotype: IgG	GenBank: NM_021962	Gene ID: 29
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29522	

Detection Antibody Information

Catalog Number: 84011-3-PBS	Clone ID: 241121G4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: active BCR-related gene
Isotype: IgG	GenBank: NM_021962	Gene ID: 29
Purification Method: Protein A purification	Immunogen Catalog Number: Ag29522	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00961-1 targets ABR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ABR Recombinant antibody, PBS Only (Capture) 84011-1-PBS (241121B4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: ABR Recombinant antibody, PBS Only (Detector) 84011-3-PBS (241121G4). 100 μ g. Concentration 1 mg/ml.

Alternative ABR matched antibody pairs: MP00961-2, MP00961-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00961-1, ABR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84011-1-PBS. Detection antibody: 84011-3-PBS. Standard: Ag29522. Range: 0.313-40 ng/mL