

For Research Use Only

# PARP14 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00966-2

## Capture Antibody Information

Catalog Number:  
84039-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241078G5

Reactivity:  
human

GenBank:  
NM\_017554

Immunogen Catalog Number:  
Ag35007

Conjugate:  
Unconjugated

Full name:  
poly (ADP-ribose) polymerase family,  
member 14

Gene ID:  
54625

## Detection Antibody Information

Catalog Number:  
84039-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241078B2

Reactivity:  
human

GenBank:  
NM\_017554

Immunogen Catalog Number:  
Ag35007

Conjugate:  
Unconjugated

Full name:  
poly (ADP-ribose) polymerase family,  
member 14

Gene ID:  
54625

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.25-80 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP00966-2 targets PARP14 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PARP14 Recombinant antibody, PBS Only (Capture) 84039-3-PBS (241078G5). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: PARP14 Recombinant antibody, PBS Only (Detector) 84039-1-PBS (241078B2). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00966-2, PARP14 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84039-3-PBS. Detection antibody: 84039-1-PBS. Standard: Ag35007. Range: 1.25-80 ng/mL