

# KIF5A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00967-1

## Capture Antibody Information

<b>Catalog Number:</b> 84013-4-PBS	<b>Clone ID:</b> 241124H1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human, mouse, rat	<b>Full name:</b> kinesin family member 5A
<b>Isotype:</b> IgG	<b>GenBank:</b> BC146670	<b>Gene ID:</b> 3798
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15518	

## Detection Antibody Information

<b>Catalog Number:</b> 84013-3-PBS	<b>Clone ID:</b> 241124F12	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> kinesin family member 5A
<b>Isotype:</b> IgG	<b>GenBank:</b> BC146670	<b>Gene ID:</b> 3798
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15518	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.313-40 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00967-1 targets KIF5A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIF5A Recombinant antibody, PBS Only (Capture) 84013-4-PBS (241124H1). 100 µg. Concentration 1 mg/ml.

Detection antibody: KIF5A Recombinant antibody, PBS Only (Capture/Detector) 84013-3-PBS (241124F12). 100 µg. Concentration 1 mg/ml.

Alternative KIF5A matched antibody pairs: MP00967-2, MP00967-3, MP00967-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00967-1, KIF5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84013-4-PBS. Detection antibody: 84013-3-PBS. Standard: Ag15518. Range: 0.313-40 ng/mL