

# PRC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00970-2

## Capture Antibody Information

**Catalog Number:**  
84041-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241216H2  
**Reactivity:**  
human  
**GenBank:**  
BC003138  
**Immunogen Catalog Number:**  
Ag8028

**Conjugate:**  
Unconjugated  
**Full name:**  
protein regulator of cytokinesis 1  
**Gene ID:**  
9055

## Detection Antibody Information

**Catalog Number:**  
84041-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241216H5  
**Reactivity:**  
human  
**GenBank:**  
BC003138  
**Immunogen Catalog Number:**  
Ag8028

**Conjugate:**  
Unconjugated  
**Full name:**  
protein regulator of cytokinesis 1  
**Gene ID:**  
9055

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00970-2 targets PRC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PRC1 Recombinant antibody, PBS Only (Capture) 84041-2-PBS (241216H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: PRC1 Recombinant antibody, PBS Only (Detector) 84041-3-PBS (241216H5). 100 µg. Concentration 1 mg/ml.

Alternative PRC1 matched antibody pairs: MP00970-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard: Ag8028. Range: 0.625-80 ng/mL

Cytometric bead array standard curve of MP00970-2, PRC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84041-2-PBS. Detection antibody: 84041-3-PBS. Standard:null. Range: null.