

FAP Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00971-1

Capture Antibody Information

Catalog Number:
84018-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241129E7

Reactivity:
human

GenBank:
NM_004460.4

Conjugate:
Unconjugated

Full name:
fibroblast activation protein, alpha

Gene ID:
2191

Detection Antibody Information

Catalog Number:
84018-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241129F6

Reactivity:
human

GenBank:
NM_004460.4

Conjugate:
Unconjugated

Full name:
fibroblast activation protein, alpha

Gene ID:
2191

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00971-1 targets FAP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FAP Recombinant antibody, PBS Only (Capture) 84018-2-PBS (241129E7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: FAP Recombinant antibody, PBS Only (Capture/Detector) 84018-3-PBS (241129F6). 100 μ g. Concentration 1 mg/mL.

Alternative FAP matched antibody pairs: MP00971-2, MP00971-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00971-1, FAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84018-2-PBS. Detection antibody: 84018-3-PBS. Standard: Eg0821. Range: 0.78-100 ng/mL