

For Research Use Only

SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00976-3

Capture Antibody Information

Catalog Number:
84060-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241215F8

Reactivity:
human

GenBank:
BC106001

Immunogen Catalog Number:
Ag13476

Conjugate:
Unconjugated

Full name:
solute carrier family 9
(sodium/hydrogen exchanger),
member 3 regulator 2

Gene ID:
9351

Detection Antibody Information

Catalog Number:
84060-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241215B6

Reactivity:
human

GenBank:
BC106001

Immunogen Catalog Number:
Ag13476

Conjugate:
Unconjugated

Full name:
solute carrier family 9
(sodium/hydrogen exchanger),
member 3 regulator 2

Gene ID:
9351

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00976-3 targets SLC9A3R2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC9A3R2 Recombinant antibody, PBS Only (Capture) 84060-4-PBS (241215F8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SLC9A3R2 Recombinant antibody, PBS Only (Detector) 84060-2-PBS (241215B6). 100 μ g. Concentration 1 mg/ml.

Alternative SLC9A3R2 matched antibody pairs: MP00976-1, MP00976-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00976-3, SLC9A3R2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84060-4-PBS. Detection antibody: 84060-2-PBS. Standard: Ag13476. Range: 1.25-40 ng/mL.