

For Research Use Only

MIRO1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01000-3

Capture Antibody Information

Catalog Number:
84055-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241157F5

Reactivity:
human

GenBank:
BC060781

Immunogen Catalog Number:
Ag15605

Conjugate:
Unconjugated

Full name:
ras homolog gene family, member T1

Gene ID:
55288

Detection Antibody Information

Catalog Number:
84055-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241157H2

Reactivity:
human

GenBank:
BC060781

Immunogen Catalog Number:
Ag15605

Conjugate:
Unconjugated

Full name:
ras homolog gene family, member T1

Gene ID:
55288

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01000-3 targets MIRO1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MIRO1 Recombinant antibody, PBS Only (Capture/Detector) 84055-2-PBS (241157F5). 100 μ g. Concentration 1 mg/ml.

Detection antibody: MIRO1 Recombinant antibody, PBS Only (Detector) 84055-4-PBS (241157H2). 100 μ g. Concentration 1 mg/ml.

Alternative MIRO1 matched antibody pairs: MP01000-1, MP01000-2, MP01000-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01000-3, MIRO1/RHOT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84055-2-PBS. Detection antibody: 84055-4-PBS. Standard: Ag15605. Range: 1.25-40 ng/mL.