

For Research Use Only

TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01006-3

Capture Antibody Information

Catalog Number:
84089-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241154B11
Reactivity:
human
GenBank:
BC048769
Immunogen Catalog Number:
Ag5792

Conjugate:
Unconjugated
Full name:
transcription factor 7 (T-cell specific, HMG-box)
Gene ID:
6932

Detection Antibody Information

Catalog Number:
84089-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241154C6
Reactivity:
human
GenBank:
BC048769
Immunogen Catalog Number:
Ag5792

Conjugate:
Unconjugated
Full name:
transcription factor 7 (T-cell specific, HMG-box)
Gene ID:
6932

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01006-3 targets TCF1/TCF7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: TCF1/TCF7 Recombinant antibody, PBS Only (Capture/Detector) 84089-2-PBS (241154B11). 100 μ g. Concentration 1 mg/ml.

Detection antibody: TCF1/TCF7 Recombinant antibody, PBS Only (Detector) 84089-3-PBS (241154C6). 100 μ g. Concentration 1 mg/ml.

Alternative TCF1/TCF7 matched antibody pairs: MP01006-1, MP01006-2, MP01006-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01006-3, TCF1/TCF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84089-2-PBS. Detection antibody: 84089-3-PBS. Standard: Ag5792. Range: 1.56-20 ng/mL.