

For Research Use Only

MYO9B Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01007-1

Capture Antibody Information

Catalog Number:

84090-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241163A5

Reactivity:

human

GenBank:

BC018108

Immunogen Catalog Number:

Ag33349

Conjugate:

Unconjugated

Full name:

myosin IXB

Gene ID:

4650

Detection Antibody Information

Catalog Number:

84090-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241163E1

Reactivity:

human

GenBank:

BC018108

Immunogen Catalog Number:

Ag33349

Conjugate:

Unconjugated

Full name:

myosin IXB

Gene ID:

4650

Applications

Tested Applications:

Cytometric bead array

Range:

1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01007-1 targets MYO9B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MYO9B Recombinant antibody, PBS Only (Capture) 84090-1-PBS (241163A5). 100 µg. Concentration 1 mg/ml.

Detection antibody: MYO9B Recombinant antibody, PBS Only (Detector) 84090-3-PBS (241163E1). 100 µg. Concentration 1 mg/ml.

Alternative MYO9B matched antibody pairs: MP01007-2, MP01007-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01007-1, MYO9B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84090-1-PBS. Detection antibody: 84090-3-PBS. Standard: Ag33349. Range: 1.563-200 ng/mL.