

For Research Use Only

MYO9B Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01007-3

Capture Antibody Information

Catalog Number:

84090-5-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241163B6

Reactivity:

human

GenBank:

BC018108

Immunogen Catalog Number:

Ag33349

Conjugate:

Unconjugated

Full name:

myosin IXB

Gene ID:

4650

Detection Antibody Information

Catalog Number:

84090-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

241163E1

Reactivity:

human

GenBank:

BC018108

Immunogen Catalog Number:

Ag33349

Conjugate:

Unconjugated

Full name:

myosin IXB

Gene ID:

4650

Applications

Tested Applications:

Sandwich ELISA

Range:

1.56-50 ng/mL (Sandwich ELISA)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01007-3 targets MYO9B in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: MYO9B Recombinant antibody, PBS Only (Capture) 84090-5-PBS (241163B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: MYO9B Recombinant antibody, PBS Only (Detector) 84090-3-PBS (241163E1). 100 µg. Concentration 1 mg/mL.

Alternative MYO9B matched antibody pairs: MP01007-1, MP01007-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology.

Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01007-3, Human MYO9B Recombinant Matched Antibody Pair - PBS only. 84090-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag33349. 84090-3-PBS was HRP conjugated as the detection antibody. Range: 1.56-50 ng/mL.



The mean MYO9B concentration was determined to be 38.37 ng/mL in K-562 cell extract based on a 2.10 mg/mL extract load.