

TRPV4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01010-2

Capture Antibody Information

Catalog Number:
84091-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241143D9
Reactivity:
human
GenBank:
NM_001177428.1
Immunogen Catalog Number:
Ag27253

Conjugate:
Unconjugated
Full name:
transient receptor potential cation channel, subfamily V, member 4
Gene ID:
59341

Detection Antibody Information

Catalog Number:
84091-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241143D2
Reactivity:
human
GenBank:
NM_001177428.1
Immunogen Catalog Number:
Ag27253

Conjugate:
Unconjugated
Full name:
transient receptor potential cation channel, subfamily V, member 4
Gene ID:
59341

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01010-2 targets TRPV4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRPV4 Recombinant antibody, PBS Only (Capture) 84091-3-PBS (241143D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRPV4 Recombinant antibody, PBS Only (Capture/Detector) 84091-2-PBS (241143D2). 100 µg. Concentration 1 mg/ml.

Alternative TRPV4 matched antibody pairs: MP01010-1, MP01010-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01010-2, TRPV4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84091-3-PBS. Detection antibody: 84091-2-PBS. Standard: Ag27253. Range: 0.781-100 ng/mL.