

For Research Use Only

Mouse Fas/CD95 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01025-2

Capture Antibody Information

Catalog Number:
84110-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241229C11
Reactivity:
mouse
GenBank:
NM_007987.2

Conjugate:
Unconjugated
Full name:
Fas (TNF receptor superfamily member 6)
Gene ID:
14102

Detection Antibody Information

Catalog Number:
84110-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
241229D3
Reactivity:
mouse
GenBank:
NM_007987.2

Conjugate:
Unconjugated
Full name:
Fas (TNF receptor superfamily member 6)
Gene ID:
14102

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01025-2 targets Fas/CD95 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Fas/CD95 Recombinant antibody, PBS Only (Capture/Detector) 84110-1-PBS (241229C11). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse Fas/CD95 Recombinant antibody, PBS Only (Detector) 84110-2-PBS (241229D3). 100 µg. Concentration 1 mg/ml.

Alternative Fas/CD95 matched antibody pairs: MP00430-1, MP00430-2, MP00430-3, MP01025-1, MP01025-3, MP50205-1, MP50205-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01025-2, MOUSE FAS/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84110-1-PBS. Detection antibody: 84110-2-PBS. Standard: Eg1375. Range: 0.625-80 ng/mL