

For Research Use Only

# Mouse Fas/CD95 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01025-3

## Capture Antibody Information

**Catalog Number:**  
84110-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241229D7  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_007987.2

**Conjugate:**  
Unconjugated  
**Full name:**  
Fas (TNF receptor superfamily member 6)  
**Gene ID:**  
14102

## Detection Antibody Information

**Catalog Number:**  
84110-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
241229D3  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_007987.2

**Conjugate:**  
Unconjugated  
**Full name:**  
Fas (TNF receptor superfamily member 6)  
**Gene ID:**  
14102

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01025-3 targets Fas/CD95 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Fas/CD95 Recombinant antibody, PBS Only (Capture) 84110-3-PBS (241229D7). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse Fas/CD95 Recombinant antibody, PBS Only (Detector) 84110-2-PBS (241229D3). 100 µg. Concentration 1 mg/ml.

Alternative Fas/CD95 matched antibody pairs: MP00430-1, MP00430-2, MP00430-3, MP01025-1, MP01025-2, MP50205-1, MP50205-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01025-3, MOUSE FAS/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84110-3-PBS. Detection antibody: 84110-2-PBS. Standard: Eg1375. Range: 1.25-80 ng/mL