

For Research Use Only

# CFB Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01031-1

## Capture Antibody Information

<b>Catalog Number:</b> 84101-3-PBS	<b>Clone ID:</b> 241141G8	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> complement factor B
<b>Isotype:</b> IgG	<b>GenBank:</b> BC007990	<b>Gene ID:</b> 629
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag0226	

## Detection Antibody Information

<b>Catalog Number:</b> 84101-1-PBS	<b>Clone ID:</b> 241141A3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> complement factor B
<b>Isotype:</b> IgG	<b>GenBank:</b> BC007990	<b>Gene ID:</b> 629
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag0226	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.781-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP01031-1 targets CFB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CFB Recombinant antibody, PBS Only (Capture) 84101-3-PBS (241141G8). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: CFB Recombinant antibody, PBS Only (Capture/Detector) 84101-1-PBS (241141A3). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative CFB matched antibody pairs: MP01031-2, MP01031-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01031-1, Complement factor B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84101-3-PBS. Detection antibody: 84101-1-PBS. Standard: Ag0226. Range: 0.781-100 ng/mL.