

PELI2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01032-1

Capture Antibody Information

Catalog Number:
84120-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241228D7

Reactivity:
human

GenBank:
BC009476

Immunogen Catalog Number:
Ag8940

Conjugate:
Unconjugated

Full name:
pellino homolog 2 (Drosophila)

Gene ID:
57161

Detection Antibody Information

Catalog Number:
84120-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241228A3

Reactivity:
human

GenBank:
BC009476

Immunogen Catalog Number:
Ag8940

Conjugate:
Unconjugated

Full name:
pellino homolog 2 (Drosophila)

Gene ID:
57161

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01032-1 targets PELI2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PELI2 Recombinant antibody, PBS Only (Capture) 84120-3-PBS (241228D7). 100 µg. Concentration 1 mg/ml.

Detection antibody: PELI2 Recombinant antibody, PBS Only (Capture/Detector) 84120-1-PBS (241228A3). 100 µg. Concentration 1 mg/ml.

Alternative PELI2 matched antibody pairs: MP01032-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01032-1, PELI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84120-3-PBS. Detection antibody: 84120-1-PBS. Standard: Ag8940. Range: 0.781-100 ng/mL.