

For Research Use Only

# ZNF350 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01034-2

## Capture Antibody Information

Catalog Number:  
84121-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241225B10

Reactivity:  
human

GenBank:  
BC009921

Immunogen Catalog Number:  
Ag23982

Conjugate:  
Unconjugated

Full name:  
zinc finger protein 350

Gene ID:  
59348

## Detection Antibody Information

Catalog Number:  
84121-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241225A9

Reactivity:  
human

GenBank:  
BC009921

Immunogen Catalog Number:  
Ag23982

Conjugate:  
Unconjugated

Full name:  
zinc finger protein 350

Gene ID:  
59348

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01034-2 targets ZNF350 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZNF350 Recombinant antibody, PBS Only (Capture) 84121-3-PBS (241225B10). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: ZNF350 Recombinant antibody, PBS Only (Detector) 84121-1-PBS (241225A9). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative ZNF350 matched antibody pairs: MP01034-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01034-2, ZNF350 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84121-3-PBS. Detection antibody: 84121-1-PBS. Standard: Ag23982. Range: 0.781-100 ng/mL.