

DPP10 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01035-2

Capture Antibody Information

Catalog Number:
84127-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241446A1

Reactivity:
human

GenBank:
BC030832

Immunogen Catalog Number:
Ag3231

Conjugate:
Unconjugated

Full name:
dipeptidyl-peptidase 10

Gene ID:
57628

Detection Antibody Information

Catalog Number:
84127-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
241446H3

Reactivity:
human

GenBank:
BC030832

Immunogen Catalog Number:
Ag3231

Conjugate:
Unconjugated

Full name:
dipeptidyl-peptidase 10

Gene ID:
57628

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01035-2 targets DPP10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DPP10 Recombinant antibody, PBS Only (Capture) 84127-1-PBS (241446A1). 100 µg. Concentration 1 mg/ml.

Detection antibody: DPP10 Recombinant antibody, PBS Only (Detector) 84127-4-PBS (241446H3). 100 µg. Concentration 1 mg/ml.

Alternative DPP10 matched antibody pairs: MP01035-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01035-2, DPP10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84127-1-PBS. Detection antibody: 84127-4-PBS. Standard: Ag3231. Range: 0.625-80 ng/mL