

For Research Use Only

# TMIGD1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP01037-2

## Capture Antibody Information

Catalog Number:  
84130-1-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241399E8

Reactivity:  
human

GenBank:  
BC137201

Immunogen Catalog Number:  
Ag25691

Conjugate:  
Unconjugated

Full name:  
transmembrane and immunoglobulin domain containing 1

Gene ID:  
388364

## Detection Antibody Information

Catalog Number:  
84130-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
241399G9

Reactivity:  
human

GenBank:  
BC137201

Immunogen Catalog Number:  
Ag25691

Conjugate:  
Unconjugated

Full name:  
transmembrane and immunoglobulin domain containing 1

Gene ID:  
388364

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP01037-2 targets TMIGD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TMIGD1 Recombinant antibody, PBS Only (Capture) 84130-1-PBS (241399E8). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: TMIGD1 Recombinant antibody, PBS Only (Capture/Detector) 84130-4-PBS (241399G9). 100  $\mu$ g. Concentration 1 mg/ml.

Alternative TMIGD1 matched antibody pairs: MP01037-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01037-2, TMI GD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84130-1-PBS. Detection antibody: 84130-4-PBS. Standard: Ag25691. Range: 1.563-200 ng/mL