

EIF3F Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP01038-4

Capture Antibody Information

Catalog Number:
84131-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240894C4
Reactivity:
human
GenBank:
BC000490
Immunogen Catalog Number:
Ag34849

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 3, subunit F
Gene ID:
8665

Detection Antibody Information

Catalog Number:
84131-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240894D2
Reactivity:
human
GenBank:
BC000490
Immunogen Catalog Number:
Ag34849

Conjugate:
Unconjugated
Full name:
eukaryotic translation initiation factor 3, subunit F
Gene ID:
8665

Applications

Tested Applications:
Sandwich ELISA

Range:
0.469-30 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01038-4 targets EIF3F in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: EIF3F Recombinant antibody, PBS Only (Capture/Detector) 84131-1-PBS (240894C4). 100 µg. Concentration 1 mg/ml.

Detection antibody: EIF3F Recombinant antibody, PBS Only (Capture/Detector) 84131-2-PBS (240894D2). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP01038-4, Human EIF3F Recombinant Matched Antibody Pair - PBS only. 84131-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34849. 84131-2-PBS was HRP conjugated as the detection antibody. Range: 0.469-30 ng/mL.



The mean EIF3F concentration was determined to be 225.37 ng/mL in Jurkat cell extract based on a 1.60 mg/mL extract load.